


JIS Type K 4.5kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Standard Blue Connector



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0677-001
General Description	4.5kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	K JIS
Magnet	30 x 20 x 20mm (L x H x W). Max Pull: 4.5kg
Cable Length	2m
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Standard Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)