



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

JIS Type K 9kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature Plug - 1.5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0678-001
General Description	Rapid, reliable temperature sensing on ferrous surfaces
Sensor Type	K JIS
Thermocouple Junction	Grounded at Tip
Magnet	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
Cable Length	1.5m
Cable Type	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
Termination Type	Miniature Plug
Max. Temperature	250°C
Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)