




# LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

## JIS Type K 9kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Standard Blue 4.5m



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

## Specifications

**Specifications**

<b>Product Code</b>	XE-0680-001
<b>General Description</b>	Rapid, reliable temperature sensing on ferrous surfaces
<b>Sensor Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded at tip
<b>Magnet</b>	25 x 25 x 40mm (HxWxL) with powerful 9kg pull
<b>Cable Length</b>	1.5m
<b>Cable Type</b>	PFA Teflon® insulated with stainless steel over-braid, fitted plug termination
<b>Termination Type</b>	Standard Plug
<b>Max. Temperature</b>	250°C
<b>Min. Temperature</b>	-50°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)