


JIS Type K 11.8kg Pull Power (Horseshoe) Magnet Thermocouple, PFA Insulated Cable with Stainless Steel Over-Braid Terminating in Miniature or Standard Plug



[Click here to view product](#)

Product Code	XE-0682
General Description	11.8kg Pull Power (Horseshoe) Magnet Thermocouple
Sensor Type	Magnet Thermocouple
Thermocouple Type	K
Magnet	45 x 30 x 30mm (L x H x W) - Max pull: 11.8 KG
Cable Length	1.5m or 3.0m
Cable Type	PFA Insulated Cable with Stainless Steel Over-Braid
Core / Strands	7/0.2mm
Cable Termination	Miniature or Standard Plug
Max. Temperature	+250°C

Min. Temperature	-50°C
Standards Met	JIS (Japanese Industrial Standards)

