



LABFACILITY

TEMPERATURE & PROCESS TECHNOLOGY

Leaf Thermocouple, Glassfibre insulated Cable with Stainless Steel Over-Braid - Type K JIS with Miniature Plug Termination - 5m length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Specifications**Specifications**

Product Code	XE-0697-001
General Description	Thermocouple mounted on a stainless steel shim
Thermocouple Type	K JIS
Thermocouple Junction	Junction grounded at tip for fast response
Cable Length	5M
Cable Type	Glassfibre insulated lead with stainless steel overbraid
Core / Strands	7/0.2mm
Termination Type	Miniature Plug (with cable clamps)
Max. Temperature	+350°C
Min. Temperature	-60°C
Standards Met	JIS (Japanese Industrial Standards)
Dimensions	25x13mm shim (0.2mm thick)