

## Nozzle Thermocouple, Glassfibre insulated Cable with Stainless Steel Over-Braid - Type K JIS with Standard Plug Termination - 2m Length



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Specifications****Specifications**

<b>Product Code</b>	XE-0891-001
<b>General Description</b>	Thermocouple housed in a stainless steel threaded bolt
<b>Thermocouple Type</b>	K JIS
<b>Thermocouple Junction</b>	Grounded hot junction for fast response, M6x1mm thread
<b>Cable Length</b>	2M
<b>Cable Type</b>	Glassfibre insulated lead with stainless steel overbraid
<b>Core / Strands</b>	7/0.2mm
<b>Termination Type</b>	Standard Plug (with cable clamps)
<b>Max. Temperature</b>	+350°C
<b>Min. Temperature</b>	-60°C
<b>Standards Met</b>	JIS (Japanese Industrial Standards)