



## PTFE Olives



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

PTFE Olives. Available to suit probe diameters ranging from 1.5mm to 8.0mm

## Specifications

<b>General Description</b>	Olives
<b>Material</b>	PTFE
<b>Probe Diameter</b>	Suitable for probe diameters between 1.5mm & 8.0mm