

Pt100 4 wire class B Resistance Thermometer with 1/2" BSP Process Connection, KNE Head



[Click here to view product](#)

Product Code	XE-3649
Sensor Type	Pt100 (100 Ohms @ 0°C) to IEC 751, Class B, 4 wire
Probe Diameter	6mm
Probe Length	150mm, 250mm, 500mm
Sheath Material	316 stainless steel
Head Type	KNE
Termination Type	4 wire connection block, M20 x 1.5mm cable entry (gland included)
Process Connection options	1/2" BSP parallel
Max. Temperature	+450°C (connection head @ 170°C)
Min. Temperature	-50°C

IP Rating	IP67
Accuracy	Class B
Number of Wires	4