

ANSI Fabricated Tube Thermocouples with Crimp Seal and Lead Wire Type K, J, T



[Click here to view product](#)

Product Code	XE-7936
General Description	Fabricated ANSI Thermocouple with Rigid Stainless Steel Sheath with 1 metre PFA lead
Thermocouple Type	ANSI Type K, J or T
Cable Length	1M
Cable Type	PFA
Termination Type	Stripped Wires
Sheath Length (mm)	50, 100, 150, 250
Sheath Diameter (mm)	3.0, 4.5, 6.0
Max. Temperature	250°C
Min. Temperature	-75°C

Standards Met	Class 1
Number of Wires	2

