



## Surface Mount RTD Sensor



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

A bolt-on RTD Temperature Sensor with an aluminium housing for fast thermal response and tab with 3mm mounting hole for easy installation, complete with a PFA lead and tails.

The sensor can be cemented, glued, screwed, or bolted on most flat surfaces.

It has a Pt100 Class A sensor as standard with lead lengths to suit, as well as non-standard connectors options, such as

[Lemo plug](#)

to suit your application.

Can also be adapted to use as a Pt1000 sensor type -  
for your custom requirements.

[contact our sales team](#)

### Specifications

## Specifications

<b>General Description</b>	Stripped 75mm tails as standard
<b>Sensor Type</b>	Pt100
<b>Fitting</b>	3mm mounting hole
<b>Cable Type</b>	PFA
<b>Termination Type</b>	Bare Tails or Miniature 3-Pin Plug Options
<b>Max. Temperature</b>	230°C
<b>Min. Temperature</b>	-50°C
<b>Tolerance</b>	Class A
<b>Length</b>	1M, 2M or 5M Lead
<b>Number of Wires</b>	4 wire or 3 wire with connector