


## WT9065 Temperature and Humidity Meter



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

The WT9065 Temperature & Humidity Meter is small and portable, with safe and reliable performance stability, a low power consumption and high precision.

The device is designed for measurement of manufacturing, laboratory, logistics, food and hygiene, air-conditioning engineering and environment protection.

Data Hold

Maximum / Minimum / Difference value

Manual / auto shutoff

LCD backlight

Power outage memory function

Temperature unit selection

### Specifications

## Specifications

<b>Product Code</b>	XE-4009-001
<b>Model</b>	WT9065
<b>Type</b>	Temperature & Humidity Meter
<b>Accuracy</b>	+/- 1°C / +/- 4% RH
<b>Resolution</b>	0.1°C (0.1°F) / 0.1% RH
<b>Battery</b>	3 x 1.5v AAA
<b>Dimensions</b>	144mm x 54mm x 32mm
<b>Humidity range</b>	5-98% RH
<b>Temperature Range</b>	-10°C to 50°C (14°F to 122°F)
<b>Weight</b>	106.1g