

Miniature Quick Wire Thermocouple Connector Socket AM-R/S-FQ Type R/S ANSI



Miniature in-line Socket - Quick Wire
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Miniature Quick Wire Connector Thermocouple Socket AM-R/S-FQ Type R/S ANSI

Easy jab-in connection, just push in wire & tighten screw. Contacts are polarised to prevent incorrect termination.

Specifications**Specifications**

Product Code	XE-1322-001
General Description	Easy jab-in connection, just push in wire & tighten screw.
Thermocouple Type	R/S ANSI
Connector Type	Quick Wire Socket
Connector Size	Miniature
Colour	Green
Max. Temperature	220°C
+Positive Contact	Copper
-Negative Contact	Copper Alloy
Standards Met	BS EN 50212:1997. RoHS