


## Standard Rectangular Fascia Thermocouple Connector Socket AS-N-FF Type N ANSI



**Standard Fascia Socket (with nylon mounting clip)**  
(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Standard Rectangular Thermocouple Connector Fascia Socket AS-N-FF Type N ANSI**

Nylon clip for quick, simple fixing to panel. All fascia mounting sockets are polarised by means of unequally sized pins

**Specifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-1577-001   |
| <b>General Description</b> | Accepts standard round pin plug connector. Nylon clip for quick, simple fixing to panel |
| <b>Thermocouple Type</b>   | N ANSI  |
| <b>Connector Type</b>      | Fascia Mounting Socket  |
| <b>Connector Size</b>      | Standard  |
| <b>Colour</b>              | Orange  |
| <b>Max. Temperature</b>    | 120°C   |
| <b>+Positive Contact</b>   | Nickel / Chromium / Silicon   |
| <b>-Negative Contact</b>   | Nickel Silicon  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS  |