


## Miniature Quick Wire Thermocouple Connector Socket JM-B-FQ Type B JIS



**Miniature in-line Socket - Quick Wire**  
(dimensions in mm)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Quick Wire Connector Thermocouple Socket JM-B-FQ Type B JIS**

Easy jab-in connection, just push in wire & tighten screw. Contacts are polarised to prevent incorrect termination.

**Specifications****Specifications**

|                            |  |
|----------------------------|--|
| <b>Product Code</b>        | XE-9223-001  |
| <b>General Description</b> | Easy jab-in connection, just push in wire & tighten screw. |
| <b>Thermocouple Type</b>   | B JIS  |
| <b>Connector Type</b>      | Quick Wire Socket  |
| <b>Connector Size</b>      | Miniature  |
| <b>Colour</b>              | Grey   |
| <b>Max. Temperature</b>    | 220°C  |
| <b>+Positive Contact</b>   | Copper   |
| <b>-Negative Contact</b>   | Copper   |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                                     |