


# Miniature Quick Wire Panel Mount Thermocouple Connector with Stainless Steel Bracket JM-K-SSPFQ Type K JIS



**Miniature Socket -**

**Quick Wire with stainless steel mounting bracket**

(dimensions in mm)



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

**Miniature Quick Wire Thermocouple Connector Panel Mount with Stainless Steel Bracket JM-K-SSPFQ Type K JIS**

Miniature size sockets supplied with stainless steel brackets, allows quick, reliable, rugged fitting into panel

**Specifications****Specifications**

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | XE-9159-001   |
| <b>General Description</b> | Contacts are polarized to ensure correct connection |
| <b>Thermocouple Type</b>   | K JIS   |
| <b>Connector Type</b>      | Quick Wire Socket with SS Bracket                   |
| <b>Connector Size</b>      | Miniature   |
| <b>Colour</b>              | Blue  |
| <b>Max. Temperature</b>    | 220°C   |
| <b>+Positive Contact</b>   | Nickel Chromium                                     |
| <b>-Negative Contact</b>   | Nickel Alloy  |
| <b>Standards Met</b>       | BS EN 50212:1997. RoHS                              |