


Barrier Terminal Block Thermocouple Connector JS-K-BTS Type K JIS



Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Barrier Terminal Block Thermocouple Connector JS-K-BTS Type K JIS

Supplied in single, interlocking pairs. Connections in thermocouple material

Specifications**Specifications**

Product Code	XE-1960-001
General Description	Supplied in single, interlocking pairs. Connections in thermocouple material
Thermocouple Type	K JIS
Connector Type	Barrier Terminal Blocks
Colour	Blue
Dimensions	See Schematic Drawing (Product Literature Tab)
Max. Temperature	220°C
+Positive Contact	Nickel Chromium
-Negative Contact	Nickel Alloy
Standards Met	BS EN 50212:1997. RoHS