




**LABFACILITY**

TEMPERATURE & PROCESS TECHNOLOGY

## Type K Retractable Curly Thermocouple Lead (IEC)




XE-7030-001

**TYPE-K (IEC)**



XE-7030-001

**TYPE-K (IEC)**



XE-7030-001

**TYPE-K (IEC)**

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Available in type K

Colour coded to IEC

2 Metres in length when fully extended

Terminated with a fitted miniature plug

Polyurethane outer

## Specifications

<b>Product Code</b>	XE-7030-001
<b>Cable Length</b>	2M fully extended
<b>Cable Type</b>	Polyurethane outer
<b>Connector Type</b>	Thermocouple Plug
<b>Connector Size</b>	Miniature
<b>Colour</b>	Green
<b>Standards Met</b>	Colour code (cores & jacket) to IEC-584-3
<b>Thermocouple Type</b>	K
<b>Termination Type</b>	Miniature Thermocouple Plug