


Type T Retractable Curly Thermocouple Lead (IEC)



XE-7033-001
TYPE-T
(IEC)



XE-7033-001
TYPE-T (IEC)



XE-7033-001
TYPE-T (IEC)

Labfacility are the UK's leading manufacturer of Temperature Sensors, Thermocouple Connectors and associated Temperature Instrumentation and stockings of Thermocouple Cables. The Company has been trading since 1971 and is ISO9001 accredited.

Available in type T

Colour coded to IEC

2 Metres in length when fully extended

Terminated with a fitted miniature plug

Polyurethane outer

Specifications

Product Code	XE-7033-001
Cable Length	2M fully extended
Cable Type	Polyurethane outer
Connector Type	Thermocouple Plug
Connector Size	Miniature
Colour	Brown
Standards Met	Colour code (cores & jacket) to IEC-584-3
Thermocouple Type	T
Termination Type	Miniature Thermocouple Plug